



# PRIORITY ONE

Independent Living and Mobility Aids



[VIEW OUR RANGE](#)

**icare**  
Medical Group

## Icare Knee Cushion



## SPECIFICATIONS

- » ActiveX™ core
- » Zip off washable cover

Measurements:

Width: 20cm

Length: 25cm

Depth: 15-7cmv

### CODES & SIZES



**Icare Knee Cushion**

Code: IC173

### Pressure Redistribution

The overlay features the ActiveX™ Lite material is the main comfort layer in this mattress. ActiveX™ is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX™ will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely.

### Ventilated Base Fabric

Vant fabric is a three-dimensional polypropylene material that has high air permeability and is anti-bacterial. Using this in the base greatly increases air movement.

#### OUR STORES

11 Formby Road  
Devonport TAS 7310

1 Edwardes Street  
Burnie, TAS 7320


#### OPENING HOURS


Monday – Friday 9:00am to 4:00pm


Other times by appointment

Closed Weekends & Public Holidays

#### CONTACT

 03 6424 5414 (DEVONPORT)

 03 6431 7752 (BURNIE)

 [priorityonemobility@bigpond.com](mailto:priorityonemobility@bigpond.com)